

MARSDEN

M-605 Large Hoist Weighing Attachment



- Class III and MDD Approved
- Accurate to 100g
- Larger than the M-600 to suit larger hoists
- Attaches to coat hanger style spreader bars
- Full numeric keypad to calculate BMI
- Optional Wi-Fi or Bluetooth available

Marsden M-605:

Large Hoist Weighing Attachment

- Class III and MDD Approved
- Accurate to 100g
- Attaches to coat hanger style weighing attachments
- Lightweight and easily portable
- Optional carry case available
- Powered by rechargeable battery pack
- 55 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Hold and Tare functions to enable accurate and easy weighing
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 8 year warranty



Ideal for care homes

The Marsden M-605 is a hoist weighing attachment which is perfect to weigh immobile patients in care homes, hospitals and private residences.

Simple hook the M-605 to your existing hoist with a coat-hanger style attachments and you are ready to begin gaining accurate readings for your patients and residents. Simple and easy to use.

Optional connectivity options

The M-605 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

Easy to use

The M-605 features from our best-selling DP-3810 indicator. This device features from a full numeric keypad for easy BMI calculation and simple function keys.

When attaching the scale to the hoist, ensure the links are vertical and within a $\pm 5^\circ$ angle to obtain an accurate weight reading.

The On/Off button ensures that you have full control over the device and you can prolong the battery life. Hold enables you to gain readings for unsteady patients and Tare allows you to deduct unwanted weight.

Capacity:	200kg	Power Supply:	Mains/Rechargeable
Graduations:	100g	Battery Life:	55 hours battery life
Functions:	Hold, Portable, Tare, BMI	Approvals:	Class III and MDD Approved
Dimensions:	730mm x 175mm x 130mm	Optional Extras:	Bluetooth/Wi-Fi, Printer, Carry Case
Weight of Scale:	6kg	Warranty Period:	8 Years