

MARSDEN

M-700 3-in-1 Baby, Child and Adult Scale



- Class III Approved
- Accurate to 20g for babies and 50g for adults
- Heavy duty and robust
- Perfect for GPs and outpatient departments
- Weigh babies, children and adults with ease
- Handrails and seat for ease of use

Marsden M-700:

3-in1 Baby, Child and Adult Scale

- Class III Approved
- Accurate to 20g for babies/children, 50g for adults
- Perfect for use in a GP surgery or outpatients department
- Powered by mains or rechargeable battery
- 55 hours of continuous use
- Auto power-off function helps to save battery power
- Heavy duty and robust
- Handrails and seat to aid use
- Easy to read display
- Hold function to enable accurate and easy weighing
- BMI calculation as standard
- Covered by the Marsden 8 year warranty



Versatile weighing

The M-700 is a versatile weighing scale as it has an ability to weigh babies, toddlers, children, adults and the elderly quickly and conveniently.

The baby tray can be easily removed to ensure there is plentiful room for adults which may wish to stand and elderly patients which may be more suited seated. Robust and heavy duty to ensure longevity.

Wirelessly transfer data

Our M-700 can be sold with Bluetooth and Wi-Fi connectivity to enable seamless data transfer to EPR systems and internal databases.

Speak to us today to see how our medical devices can connect into your system and how we can work together to ensure that we can meet all of your requirements in regards to data transfer.

Easy to use

The M-700 features from our best-selling DP-3810 indicator. This device will measure a patient's weight, a full numeric keypad and a secondary display to show BMI.

Hold function allows you to accurately capture a weight, even when the patient is moving. BMI can also be easily be calculated using the full numeric keypad on the indicator.

| | | | |
|-------------------------|---|-------------------------|----------------------------|
| Capacity: | 150kg | Power Supply: | Mains/Rechargeable |
| Graduations: | Baby/Toddler: 60kg x 20g Adult: 150kg x 50g | Battery Life: | 40 hours of continuous use |
| Functions: | BMI, Portable, Optional Data Transfer | Approvals: | Class III Approved |
| Dimensions: | 700mm x 630mm x 1250mm | Optional Extras: | Bluetooth/Wi-Fi, Printer |
| Weight of Scale: | 4kg | Warranty Period: | 8 Years |