

MARSDEN

PRECISION WEIGHING SOLUTIONS

MHS-2710 Hoist Weighing Attachment

Class III Approved

Accurate to 100g

Suitable for most hoist weighers

Simple and easy to use

Lightweight and easily portable

Powered by rechargeable battery

100 hours of continuous use

Auto power-off function helps to save



Ideal for Care Homes

The Marsden MHS-2710 is the ideal solution for weighing immobile patients in care homes, hospitals, and private residences.

Designed for seamless integration, it easily attaches to your existing hoist, providing instant access to accurate weight readings with minimal effort.

Ensure comfort and precision for your patients and residents with this dependable weighing attachment.

Simple to Use

The MHS-2710 is designed with user convenience at its core. Its large, clear LCD display ensures effortless visibility of weight readings, even when mounted on a hoist.

Simplifying operation, the scale features an intuitive two-button interface, making it quick and easy to accurately measure a patient's weight.

Simple to set up

The MHS-2710 is designed with user convenience at its core. Its large, clear LCD display ensures effortless visibility of weight readings, even when mounted on a hoist.

Simplifying operation, the scale features an intuitive two-button interface, making it quick and easy to accurately measure a patient's weight..

Capacity:	300kg	Power Supply:	Rechargeable Battery
Graduations:	100g	Battery Life:	100 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Class III approved
Dimensions:	192mmx50mmx121mm	Optional Extras:	N/A
Weight:	1kg	Warranty Period:	8 years