










# VET SCALES COMPARISON TABLE

	V-20	V-22	V-25	V-100	V-110	V-150	V-180	V-250	VT-250
PRODUCT IMAGE									
PRODUCT DESCRIPTION	PREMIUM SMALL ANIMAL SCALE	SMALL VET SCALE FOR CATS AND DOGS	LOW COST SMALL VET SCALE	MEDIUM ANIMAL WEIGHING SCALE	PORTABLE, HIGH CAPACITY FOR LARGE PETS	LARGE HIGH CAPACITY VET SCALE	AFFORDABLE LARGE VET SCALE	LARGE HIGH CAPACITY STAINLESS STEEL SCALE	CONSULTING ROOM TABLE WITH SCALE
CAPACITY	20KG	20KG	25KG	100KG	150KG	150KG/300KG	180KG	250KG	250KG
GRADUATIONS	10G	5G	10G	50G	100G	50G/100G	100G	100G	100G
APPROVALS	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
WARRANTY PERIOD	1 YEAR	1 YEAR	1 YEAR	1 YEAR	1 YEAR	1 YEAR	1 YEAR	1 YEAR	1 YEAR
DIMENSIONS	620MM X 350MM X 60MM	585MM X 335MM X 115MM	630MM X 400MM X 85MM	715MM X 410MM X 40MM	900MM X 550MM X 62MM	900MM X 550MM X 80MM	900MM X 500MM X 25MM	1090MM X 505MM X 70MM	1160MM X 550MM X 830MM
WEIGHING UNITS	KG/LB	KG/LB	KG/LB	KG/LB	KG/LB	KG/LB	KG/LB	KG/LB	KG/LB
WEIGHT OF SCALE	2.5KG	2.7KG	1.8KG	8KG	5.5KG	15KG	6.2KG	18KG	35KG
POWER SUPPLY	9V BATTERY	4 X AA BATTERIES	2 X AAA BATTERIES	2 X AAA BATTERIES	2 X AAA BATTERIES	RECHARGEABLE BATTERY	4 X AA BATTERIES	6 X AA BATTERIES	4 X AA BATTERIES
BATTERY LIFE	40 HOURS	40 HOURS	40 HOURS	40 HOURS	40 HOURS	40 HOURS	40 HOURS	40 HOURS	40 HOURS
DATA TRANSFER CAPABILITY?	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
IDEAL FOR	HOME AND SURGERY	HOME AND SURGERY	HOME AND SURGERY	HOME AND SURGERY	HOME AND SURGERY	HOME AND SURGERY	HOME AND SURGERY	VET SURGERY	CONSULTING ROOM